


PRODUCT IMAGE



PRODUCT DETAILS

Total Layer : 5
Size : 60x90 cm
Weight : 60gm
Absorbent Core Material :
Super Absorbent Polymer (SAP)
Super Absorbent Polymer Content : 3gm +-1gm
Absorbent Rate : NMT 20 sec
Moisture : <_ 10%
Absorbency : Heavy Absorbency, +-500ml
Length of Filler Materials : 76cm +-1cm
Width of Filler Materials : 50cm +-1cm

APPLICATION



the people with reduced mobility



postpartum care



menstrual women



bedridden



incontinence



patient after surgery

First Layer: Soft Non Woven Fabric which covers the absorbent filler and prevent the filler from reaching the skin or clothes under normal handling

Second Layer : Tissue Paper/ Non Woven

Third Layer : Fluff Plup mixed with SAP
(Absorbent the liquid very fast and quickly)

Fourth Layer : Tissue Paper

Fifth Layer : PE Film (prevent leakage and keeps the bed or floor dry and clean)

